

Technical data sheet
Lymphocyte Separation Media
(Density 1,077 g/ml)

CAT N°: AC-AF-0017 and AC-AF-0018

Theoretical pH: 7.0 ± 0.5

Osmolality: 300 mOsm/kg ± 20

Density: 1.077 ± 0.001

Colour: colourless, clear solution

Storage conditions: Room temperature, protected from light

Shelf life: 24 months

Sterility tests:

- bacteria in aerobic and anaerobic conditions
- fungi and yeast

Endotoxin: < 10 EU/ml

Composition: Available on request

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used for in vitro diagnostic and scientific purposes. Do not use it in therapy, human or veterinary applications.

Applications:

Lymphocyte Separation Media is designed for the simple, rapid isolation of lymphocytes from whole blood that has been diluted and treated with anti-coagulant or defibrinating agent. For best results, use blood drawn less than 2 hours before. Do not use blood more than 24 hours from when it was drawn.

Uses:

- 1) Thoroughly mix the Lymphocyte Separation Media by inverting the bottle gently.
- 2) Aseptically transfer 3 ml of Lymphocyte Separation Media to a 15 ml centrifuge tube.
- 3) Mix 2 ml of defibrinated or heparinated blood with 2 ml of physiological saline (PBS w/o Ca w/o Mg) or balanced salt solution (AC-BS-0002).
- 4) Carefully layer the diluted blood over 3 ml of Lymphocyte Separation Media (room temperature) in a 15 ml centrifuge, creating a sharp blood-Lymphocyte interphase. DO NOT MIX! The quality of the separation is dependent upon a sharp interphase between the lymphocytes and the solution.
- 5) Centrifuge the tube at 400G at room temperature for 15 to 30 minutes. Centrifugation should sediment erythrocytes and polynuclear leukocytes and band mononuclear lymphocytes above the Lymphocyte Separation Media.
- 6) Aspirate the top layer of clear plasma to within 2-3 mm above the lymphocyte layer.
- 7) Aspirate the lymphocyte layer plus about half of the Lymphocyte Separation Media layer below it and transfer it to a centrifuge tube. Add an equal volume of buffered balanced salt solution to the lymphocyte layer in the centrifuge tube and centrifuge for 10 minutes at room temperature (18°C to 25°C) at a speed sufficient to sediment the cells without damage i.e., 160-260 g. Washing the cells removes Lymphocyte Separation Media and reduces the percentage of platelets.
- 8) Wash the cells again with buffered balanced salt solution (AC-BS-0002) and resuspend in the appropriate medium for your applications.

Important Remarks:

- CAUTION : the product is not for human or animal therapeutic use. Uses other than the intended use may be a violation of local law.
- Each laboratory must carry out their own testing procedures on new media according to national legislation prior to releasing them to the lab for routine in vitro applications.
- Each clinician/scientist must make an independent judgment on whether this medium is suitable for use in in vitro diagnostic applications conducted in their laboratory.
- anprotec does not guarantee the successful outcome of any diagnostic testing based solely on the use of Anprotec brand medium.