

Trypsin – EDTA (0,05%) 1X in solution
w/o Calcium w/o Magnesium w/ Phenol Red

AC-EZ-0004

Components	Quantity in g/l
EDTA Disodium Salt Dihydrate	0,22140000
Potassium Chloride	0,40000000
Sodium Chloride	8,00000000
D-Glucose Anhydrous	1,00000000
Phenol Red Sodium Salt	0,00200000
Trypsin 1:250 (Non-ACF)	0,50000000
Sodium Bicarbonate	0,58000000